



## MEASURING SCREW LENGTH BY HIT TECHNIQUE

**Dr Ashwani kumar  
Katotra**

Senior resident, post graduate department of orthopedics, GMC Jammu.

**Dr Mohinder singh  
chib\***

Senior resident, post graduate department of orthopedics, GMC Jammu.  
\*Corresponding Author

**Dr Farid H. Malik**

Lecturer, post graduate department of orthopedics, GMC Jammu.

**ABSTRACT**

It is seen that while treating fractures surgically using plating one of the challenge is to measure exact screw size to fix the fracture. Various modalities are being practiced like use of depth gauge,[1] fluoroscopic guidance and clinical experience. Our study shows that among all methods used, there is one more method that we are using since 2 years (January 2018— January 2020) has much more reliability and also no need of any external appliance to depend upon. The method named by us as HIT TECHNIQUE proved to be very much reliable in assessing exact screw length.

**KEYWORDS :****INTRODUCTION:**

Plating in orthopedic surgeries has revolutioned the branch by providing anatomical reduction and rigid fixation.[2]. Postoperative radiograph assessment mainly focuses parameters like reduction of fracture and stability of fixation. Stability is provided by bicortical purchase of screw.[3,4]. It is seen that the fear of bicortical purchase many times leads to long or short screw insertion. The error occur mostly due to misinterpretation by depth gauge. On fluoroscopic assessment if size found to be wrong then replacing by another screw make the hole wide and thus loosening of screw and affecting fracture stability.[5]

**MATERIALS AND METHODS:** The study being conducted in tertiary hospital orthopedic department. All subjects undergoing fixation using plating technique were subjected to screw length measurement using HIT TECHNIQUE. No depth gauge or fluoroscopy was used for assessing screw length. The technique was used only for assessing screw length in diaphyseal region of bone.[6,7,8] The technique can be applied in any long bone like humerus, radius, ulna, femur, tibia, metacarpal shafts. It has no role in cancellous bones and metaphyseal region of bone. Results are assessed in terms of accuracy of screw length in postoperative radiographs. Cases with pathological bones and wide osteoporosis are excluded from study.

**TECHNIQUE:**

**HIT TECHNIQUE:** It involves the measurement of screw length without the use of depth gauge or fluoroscopy. The method involves the above mentioned steps.

**Step 1;** Using the drill and sleeve, first cortex is drilled and as the second cortex is reached further drilling is stopped.

**Step 2;** If the drill bit has any markings over it read the mark on drill bit when second cortex hits. If drill bit has no markings remove the drill bit once second cortex hits.

If drill bit has markings over it then after reading mark continue to drill. Once drill bit is removed apply sleeve over drill bit and reach the marking read during hitting of second cortex. The drilling part of drill bit unexposed by drill sleeve is measured as distance X mm. (X—arbitrary used)

If drill bit has no markings over it then remove the drill bit after drilling one cortex as it hits the second cortex. Use K-wire and measure the drilled distance equal to distance of K-wire inside bone. This is equivalent to distance X.

Step 3; Now screw length can be found easily by using distance X, which is equal to X+2mm in case of small bones like metacarpals and forearm bones and X+4mm in case of bones like humerus, femur and tibia.

**DISCUSSION:**

It is seen that the technique is very reliable in calculating the screw length without use of any depth gauge and fluoroscopy. As the technique is very much sensitive with little scope of error as encountered with depth gauge due to improper fitting across bone. The technique is proved to be very effective in routine surgical workup in treating diaphyseal fractures like humerus shaft, forearm bones, metacarpals, phalanges, femur shaft and tibia diaphyseal fractures. The technique has drawback in cancellous metaphyseal zone of bone like proximal humerus (head of humerus), distal humerus, proximal ulna, distal radius, pelvic bone, trochantric area, distal femur, tibia plateau and pilon, talus and calcaneum where second hit cannot be appreciated. Surgeon precision to judge the reach of second cortex is main dependent variable. High speed drills to be avoided in procedure as using these appliances causes difficulty to appreciate second cortex.

**REFERENCES:**

1. Williams D F, Roaf R 1973. Implants in Surgery. Saunders London.
2. Anderson L D, Sisk T D, Tooms R E, Park III W 1975. Compression plate fixation in acute diaphyseal fractures of radius and ulna. J Bone Joint Surg 57 A287-297.
3. Peterson I.T. Principles of internal fixation with plates and screws. Arch Surg 64: 345-354.
4. Uhl R L 1989 the biomechanics of screws. Orthopedic Review 18; 1302-1307
5. Chao E Y S, Aro H T 1991. Biomechanics of fracture fixation. In: Mow V C, Hayes W C. Basic Orthopedic Biomechanics. Raven Press, New York pp 293—336. Müller et al., 1990.
6. Müller ME, Allgöwer M, Schneider R, et al: Manual of internal fixation: techniques recommended by the AO-ASIF group, 3rd ed. Berlin, Springer-Verlag, 1990.
7. Schatzker and Tile, 1980. Schatzker J, Tile M: The AO (ASIF) method of fracture care. Instr Course Lect 1980; 29: 41.
8. Behrens F. General Theory and Principles of External Fixation. Clin Orthop 1989; 241: 15.