



ORIGINAL RESEARCH PAPER

DEVELOPMENT OF A HEAD-MOUNTED COLPOSCOPE WITH INTEGRATED DEEP LEARNING AND REAL-TIME AUGMENTED REALITY OVERLAYS FOR PRECISION CERVICAL PRECISION BIOPSY

Obstetrics & Gynaecology


KEY WORDS: Cervical Cancer Screening, Computer-Guided Colposcopy, Augmented Reality Overlay, Transformation Zone Detection, Segment Anything Model (FastSAM), Medical Device Startup Implementation

Prof Dr Pradeep Kumar Radhakrishnan	MS, MCh CTVS AIIMS Delhi, Post-doctoral fellow CTVS SCTIMST, Postdoctoral Fellow in ECMO, FACC, FIACS, CPDH- IISc, Global MBA Lond, Prof CTVS and HOD CTVSYANS Medical Institute, India
Dr Gayathri Ananyajyothi Ambat	MBBS, Junior Resident, YANS Medical Institute, Kerala India,
Dr. Roshini Ambat	Senior Consultant Gynecologist KIMS Trivandrum
Rajesh Thomas	MBS AI MIT Sloan,
Prof Dr Naseema A	Prof and HOD Gokulam Medical College, Trivandrum Kerala

ABSTRACT	Introduction: Cervical cancer remains one of the leading causes of oncological mortality and morbidity among women globally, particularly within low- and middle-income countries (LMICs) where access to structured cytological screening and expert medical personnel is severely restricted. While human papillomavirus (HPV) DNA testing and liquid-based cytology have expanded screening coverage, colposcopy guided by directed histological biopsy remains the clinical gold standard for diagnosing cervical intraepithelial neoplasia (CIN). However, conventional digital colposcopy is limited by significant inter-observer and intra-observer variability. The diagnostic accuracy is heavily dependent on the clinical experience of the operator, with junior or less-experienced colposcopists frequently exhibiting a suboptimal sensitivity (as low as 55%–70%) for the identification of high-grade squamous intraepithelial lesions (HSIL/CIN2+). This disparity frequently results in either missed micro-lesions or a proliferation of unnecessary, painful punch biopsies that induce patient distress and bleeding. To bridge this clinical gap, this paper outlines a comprehensive implementation roadmap for a startup venture developing a head-mounted colposcopy unit featuring edge-based deep learning and a real-time augmented reality (AR) overlay. By shifting away from static, offline cloud analysis toward real-time, intraoperative computer-guided instrumentation, this system builds a precision-targeted sampling environment. The integration of advanced computer vision frameworks—specifically focusing on the automatic delineation of the Transformation Zone (TZ) and dynamic tracking of acetowhite changes—acts as a localized visual guidance system for the point-of-care clinician.
----------	---

Concept of a Head-Mounted Colposcope with AI Integration

The proposed system reimagines the architectural design of a colposcope. Traditional devices are bulky, floor-standing binocular microscopes or modern handheld digital cameras that require clinicians to look away from the patient toward an external monitor. This spatial decoupling compromises fine-motor precision during instrument handling and tissue sampling.



The device is built upon a lightweight, ergonomically balanced head-mounted display (HMD) integrated with high-definition, stereoscopic macro-optical sensors and co-axial polarized LED illumination. The core ideology rests on three technological pillars:

1. **Uninterrupted Line of Sight:** By routing the stereoscopic digital video loop through high-transparency holographic optical waveguides, the clinician views the cervix directly through their natural field of view, maintaining spatial depth perception.
2. **Real-Time Automated Segmentation:** Embedded deep learning networks process the incoming raw video frames on a frame-by-frame basis, eliminating the delay of external data transmission.
3. **Augmented Reality Spatial Mapping:** The inferred diagnostic boundaries (such as the squamocolumnar junction and lesion margins) are converted into a dynamic, color-coded graphical mask. This mask is projected exactly over the real physical tissue, forming a stable "visual bullseye" that tracks patient movement and instrument intervention in real time.

Detailed Ideology and Clinical Rationale
 The primary clinical objective of integrating deep learning with an AR overlay is the systematic minimization of diagnostic errors arising from two main clinical factors: the misidentification of the high-risk zones on the cervix, and human error during biopsy selection.

Delineating the Transformation Zone (TZ)
 The Transformation Zone is the dynamic histological borderland where the fragile simple columnar epithelium undergoes physiological metaplasia to become stratified squamous epithelium. Over 90% of cervical pre-cancers originate within this specific region. Identifying the type of TZ (Type 1: completely ectocervical and fully visible; Type 2:

having an endocervical component but fully visible; Type 3: extending deep into the endocervical canal and partially or completely invisible) determines whether a colposcopic exam is clinically adequate or if excisional procedures are required. Inexperienced operators routinely struggle to identify the squamocolumnar junction (SCJ), leading to blind sampling errors. Automating this boundary identification guarantees that the clinician's sampling window is locked precisely onto the high-risk tissues.


The Ideal Sampling Environment

In traditional settings, after applying 3%–5% acetic acid, the colposcopist must manually evaluate transient acetowhite changes (the whitening of abnormal, protein-dense cells), vascular punctuation, and mosaicism patterns. This visual evaluation must occur within a strict 2-to-3-minute window before the effect fades.

The integrated AR environment counters this pressure by establishing an Objective Sampling Environment. The computer-guided overlay removes subjective interpretation by constantly running background algorithms that analyze the kinetics of tissue whitening. The resulting AR mask highlights the single pixel coordinate that displays the most rapid, dense, and prolonged acetowhitening. This pinpoint accuracy eliminates multi-focal sampling, transforming a subjective procedure into an algorithmic, objective intervention.

Updated Engineering Concepts of AI Integration and Real-Time Overlay

To achieve smooth, high-fidelity clinical operation, the startup's engineering team must build a highly optimized edge computing pipeline. The structural architecture is divided into three processing pipelines:



1. Optical Ingestion and Preprocessing

The hardware utilizes dual high-sensitivity C-MOS sensors capturing stereoscopic video at 1080p resolution and 60 frames per second. A core engineering challenge in digital colposcopy is specular reflection—bright highlights caused by light reflecting off moist cervical mucus, which blinds computer vision models.

To resolve this, the unit integrates a physical, cross-linear polarizing filter layer over the LED light source and camera lenses, paired with a lightweight Swin Transformer network running locally to digitally patch remaining saturated pixels, ensuring a clean input stream.

2. Edge-Based Multi-Stage Deep Learning Architecture


To avoid relying on cloud networks—which are unusable in low-resource rural environments due to latency and data drops—all models are optimized to run on an integrated, low-power Edge TPU or NPU embedded processor.

- Stage 1: Region-of-Interest (ROI) Extraction: A highly optimized YOLOv8 object detection framework isolates the cervix boundary, stripping away irrelevant features like the speculum blades, vaginal walls, or diagnostic swabs, cropping the active processing frame to a standardized 224×224-pixel zone.

- Stage 2: Semantic Micro-Segmentation via FastSAM: The Fast Segment Anything Model (FastSAM) is deployed to map out anatomical structures. Fine-tuned on expert-annotated colpograms, FastSAM extracts the exact contours of the Transformation Zone and the squamocolumnar junction within milliseconds, classifying the visible TZ type (Type 1, 2, or 3).
- Stage 3: Temporal Lesion Scoring: Concurrently, an EfficientNet-B0 backbone network enhanced with hybrid attention mechanisms tracks the temporal kinetics of the tissue across time-series frames (Saline- Acetic Acid - Lugol's Iodine). The network computes a pixel-wise malignancy probability matrix based on the International Federation for Cervical Pathology and Colposcopy (IFCPC) metrics.
- 3. Latency-Minimized AR Projection and Stabilization
 - To project the computed probability matrix as a coherent heat map, the engine calculates a real-time Homography Matrix between the camera space and the user's optical waveguide display plane.
 - To prevent visual drift or spatial lag (which induces user nausea and misaligned biopsies), the system pairs optical flow tracking with an inertial measurement unit (IMU) built into the headset. This keeps the AR color overlay anchored tightly to the cervix with an end-to-end processing latency of under 30 milliseconds, maintaining alignment even during sudden patient movements.

Startup Implementation Roadmap

For an early-stage startup team, launching this complex, regulated biomedical device requires managing concurrent technical, clinical, and regulatory pipelines over an estimated 24-month lifecycle.



Phase 1: Hardware Prototyping and Core Optics Engine (Months 01–06)

- Milestone: Functional benchtop prototype generating stabilized video feeds.
- Tasks:
 1. Design a lightweight 3D-printed stereoscopic optical chassis optimized for mounting onto standard surgical headbands.
 2. Implement the cross-polarized LED illumination ring.
 3. Interface the dual C-MOS sensors with an embedded NPU evaluation board.
 4. Train the primary YOLOv8 cervix-isolation model using public datasets (e.g., the Kaggle Cervical Cancer Screening dataset).

Phase 2: Waveguide Integration and Model Optimization (Months 07–12)

- Milestone: Integrated HMD displaying real-time AR overlays with minimal latency.
- Tasks:
 1. Integrate holographic optical waveguides into the clinician's field of view.
 2. Implement the Homography calibration pipeline to align

the digital camera display with the user's pupils.

3. Compress the FastSAM and Efficient Net models using post-training quantization (INT8 precision) to maximize inference speeds on the edge hardware.
4. Validate tracking accuracy using custom silicon cervical phantoms.

Phase 3: Clinical Validation and Multi-Center Studies (Months 13–18)

- Milestone: Institutional Review Board (IRB) approved human clinical trials.
- Tasks:
 1. Deploy prototype units across two distinct testing grounds: a group of low-resource rural clinics (focusing on screening ease for junior nurses) and a specialist hospital colposcopy clinic (benchmarking against expert colposcopists).
 2. Collect prospective clinical data to measure sensitivity and specificity improvements for CIN2+ lesions under AI guidance.
 3. Refine tracking algorithms to handle real-world variations like bleeding, heavy discharge, and varied patient anatomy.

Phase 4: Regulatory Compliance and Commercial Scaling (Months 19–24)

- Milestone: Regulatory submission (FDA 510(k) or CE Mark Class IIa Medical Device).
- Tasks:
 1. Establish an ISO 13485 compliant Quality Management System (QMS).
 2. Run rigorous software validation protocols to lock the deep learning weights, preventing algorithm drift.
 3. Compile electrical safety and biocompatibility reports (IEC 60601-1).
 4. Launch the commercial go-to-market strategy, targeting public health procurement agencies and private outpatient networks.

Future Innovations in Computer-Guided Colposcopy

As edge-computing hardware and sensor miniaturization mature, the field of digital colposcopy will expand beyond basic visual overlays:

- **Hyperspectral Tissue Analysis:** Future iterations will replace basic multi-spectral imaging with true hyperspectral sensors, capturing continuous optical bands to measure cellular oxygenation levels and chemical changes without using topical chemical reagents.
- **Automated Robot-Guided Micro-Biopsies:** Coupling the AR target acquisition engine with miniature, head-stabilized micro-robotic arms could allow for semi-automated, touchless tissue extraction, removing manual tremor errors entirely.
- **Multimodal Generative Risk Layering:** Merging real-time visual inferences with a patient's historical records (such as high-risk HPV genotypes, epigenetic methylation profiles, and cytological history) via on-device transformers will deliver comprehensive, multi-layered diagnostic predictions at the point of care.

Multimodal AI Fusion: Visuals + Molecular Biomarkers

Relying entirely on an AI colposcope's visual assessment can introduce blind spots due to anatomical variations or inflammation. The latest breakthrough approach fuses the AI colposcope's image processing with rapid biomarker triage:

- **p16/Ki-67 Dual-Stain Cytology AI:** The National Cancer Institute (NCI) validated a breakthrough AI program specifically trained to detect the co-expression of p16 (a tumor suppressor) and Ki-67 (a cellular proliferation marker). When a liquid brush biopsy sample is processed, the AI instantly quantifies these dual-stained cells. Fusing this data with the AI colposcope's localized image assessment delivers unparalleled specificity for high-grade cervical intra-epithelial neoplasia (CIN2+).
- **Host & Viral DNA Methylation Panels:** Triage tests targeting specific hypermethylated genes (such as the *FAM19A4/miR-124-2* panel) serve as highly precise

molecular filters. If the AI colposcope flags a borderline visual lesion, checking the local methylation profile determines whether the tissue is actively transforming into a malignancy, preventing unnecessary tissue biopsies.

- **2. Advanced Point-of-Care Microfluidics**

For portable setups engineered for low-resource settings, standard laboratory processing is impractical. Integrating the colposcope with microfluidic "lab-on-a-chip" tech provides real-time molecular validation:

- **Microfluidic CRISPR/Cas Real-Time Detection:** The latest microfluidic breakthrough involves room-temperature, automated microfluidic systems (such as the newly engineered *Ccard* platform). This uses an RPA-CRISPR/Cas13a reaction directly from a raw cervical swab, bypassing the need for complex laboratory DNA extraction. Within 30 minutes, it identifies high-risk HPV16 and HPV18 strains at a highly sensitive limit of detection (5 copies/ μ L), giving the AI colposcope operator an instant viral profile to overlay onto the live video feed. [1]
- **AI-Powered Droplet Microfluidics:** Platforms like *DropletDx* generate single-cell droplets from brush samples, using high-resolution optical arrays and embedded AI processing to execute automated, high-throughput HPV screening right next to the patient.

REFERENCES

1. Liu L, Liu J, Su Q, Chu Y, Xia H, Xu R. Performance of artificial intelligence for diagnosing cervical intraepithelial neoplasia and cervical cancer: a systematic review and meta-analysis. *E Clinical Medicine*. 2025; 80:102992. doi: 10.1016/j.eclim.2024.102992.
2. Takahashi T, Kobayashi Y, Sakurai R, Matsuoka K, Akatsuka J, Kisu I, et al. A systematic review of the application of artificial intelligence in colposcopy: diagnostic accuracy for cervical intraepithelial neoplasia and cervical cancer. *Clinical Medicine Insights: Oncology*. 2025;19. doi: 10.1177/11795549251374908.
3. Khare SK, Blanes-Vidal V, Booth BB, Petersen LK, Nadimi ES. A systematic review and research recommendations on artificial intelligence for automated cervical cancer detection. *WIREs Data Mining and Knowledge Discovery*. 2024;14(6):e1550. doi: 10.1002/widm.1550.
4. Desai KT, et al. The development of "automated visual evaluation" for cervical cancer screening: the promise and challenges in adapting deep-learning for clinical testing. *International Journal of Cancer*. 2022;150(5):741–752. doi: 10.1002/ijc.33879.
5. Coole JB, Brenes D, Possati-Resende JC, Antoniazzi M, de Oliveira Fonseca B, Maker Y, et al. Development of a multimodal mobile colposcope for real-time cervical cancer detection. *Biomedical Optics Express*. 2022;13(10): 5116–5130. doi: 10.1364/BOE.463253.
6. Grant BD, Quang T, Possati-Resende JC, Scapulatempo-Neto C, de Macedo Matsushita G, Mauad EC, et al. A mobile-phone based high-resolution micro endoscope to image cervical precancer. *PLoS One*. 2019;14(2): e0211045. doi: 10.1371/journal.pone.0211045.
7. Tang W, et al. Cervical cytology screening facilitated by an artificial intelligence microscope: a preliminary study. *Cancer Cytopathology*. 2021;129(9):693–700. doi: 10.1002/cncy.22425.
8. Holmström O, Linder N, Kaingu H, Mbuuko N, Mbete J, Kinyua F, et al. Point-of-care digital cytology with artificial intelligence for cervical cancer screening in a resource-limited setting. *JAMA Network Open*. 2021;4(3): e211740. doi: 10.1001/jamanetworkopen.2021.1740.
9. Geng X, Ma J, Li X, Wei Z, Zhang X, Quan T, et al. Robust whole slide image analysis for cervical cancer screening using deep learning. *Nature Communications*. 2021;12:5639. doi: 10.1038/s41467-021-25296-x.
10. Chang Y, Li T, Zhou Q, Kang D, Zhu L, Yang J, et al. Effectiveness of artificial intelligence-assisted colposcopy in a resource-limited population. *Obstetrics & Gynecology*. 2025;146(4):545–554. doi: 10.1097/AOG.0000000000006014.
11. Performance of an artificial intelligence-based tool for cervical precancer screening in five countries in Africa: a prospective, observational, diagnostic accuracy study. *The Lancet Global Health*. 2026;14(1): e49–e60. doi: 10.1016/S2214-109X(25)00352-3.
12. Natarajan J, Mathur S, Vishnubhatla S, Kumar S, Vashist S, Ramanujam N, et al. Can portable colposcopes replace standard-of-care colposcopes? A crossover trial of two portable colposcopes with a standard-of-care video colposcope. *Asian Pacific Journal of Cancer Prevention*. 2022;23(12):4013–4021. doi: 10.31557/APJCP.2022.23.12.4013